IHDEDENTAL*



IMMEDIATE LOADING DENTAL IMPLANT SYSTEM


THE ADVANTAGES

OF THE ENDOSSEOUS DENTAL IMPLANT SYSTEM ROOT®

KOS **ROOT**[®] implants are used for crowns, bridges and bars. The compression screw design permits to incorporate the restoration in an immediate loading protocol (incorporation of the prosthesis within max. three days). The single-piece design reduces costs, the danger of periimplantitis and it eliminates the hazzles of screw loosening. KOS **ROOT**[®] implants are straight with prosthetic head for cementation.

The prescribed or recommended tightening torques for implantscan be found on our website:

www.implant.com/en/downloads






Wide range of sizes Enossal length 6 - 20 mm Enossal Ø 3.0 - 5.5 mm Angulation adapter 15° and 25° No-Itis® Laser Surface A smooth surface that, in contact with the bone is shaped like a rough surface



The new surface treatment for Dr. Ihde Dental AG implants is created with the latest generation of robotic tools for laser ablation. This new technology of high precision creates roughness in the implant through a mesh of hemispherical micrometric pores, with a defined, always identical size and shape and with a symmetrical distribution. The result is a more adequate topography, which provides the most suitable conditions for the osseointegration of the implant, but at the same time it is, and behaves like, a smooth surface at a micrometric (cellular) level. This means that while bone grows well on this surface, the adhesion of bacteria to the same surface is significantly reduced.



a) Abutment Ø b) Abutment height c) Neck length d) Enossal length e) Enossal Ø f) Neck Ø KOS **ROOT® 3.0 - 4.0**



3.9 mm

7.2 mm

3.7 mm

6 - 20 mm

3.0 - 5.5 mm

2.0 - 2.55 mm

Max. insertion torque 50 Ncm

Max. insertion torque 80 Ncm

KOS **ROOT**[®] implants are delivered incl. lab-set REF 462353, consisting of

Impression and castable abutment, engages in the three vertical notches

Double analogue, plastic

IA4/IAU 462111

PA X 462136







NOTE TSPA 4 is part of the standard lab-set, but CAN NOT be used for KOS **5. ROOT**[®]. **TSPA 4** can only be used on implants with a small abutment head, which can be found in the extended lhde Dental KOS[®] system application..

462029



Neck Ø 2.0 mm:



Enossal Ø		Neck Ø	Note		Price cat.
	6 mm			458100	
	8 mm			458101	
	10 mm			458102	
3.0 mm	12 mm	2 mm	bendable	458103	F
	14 mm	2	Seriodore		
	16 mm 18 mm			458105 458106	
	20 mm			458107	
	6 mm			458110	
	8 mm			458111	
	10 mm			458112	
	12 mm			458113	
3.5 mm	14 mm	2 mm	bendable		F
	14 mm			458115	
	18 mm			458116	
	20 mm			458117	
	6 mm			458120	
	8 mm			458121	
	10 mm			458122	
	12 mm			458123	
4.0 mm	14 mm	2 mm	bendable		F
	16 mm			458125	
	18 mm			458126	
	20 mm			458127	
	6 mm			458130	
	8 mm			458131	
	10 mm			458132	
	12 mm			458133	
4.5 mm	14 mm	2.35 mm		458134	F
	16 mm			458135	
	18 mm			458136	
	20 mm			458137	
	6 mm			458140	
	8 mm			458141	
5.0 mm	10 mm	2.35 mm		458142	F
	12 mm			458143	
	14 mm			458144	
	6 mm			458150	
	8 mm			458151	
5.5 mm	10 mm	2.55 mm		458152	F
	12 mm			458153	
	14 mm			458154	

ADDITIONAL IMPRESSION TAKING AND LABORATORY ACCESSORIES



Description	Unit	Material	Code	REF	Price cat.
Impression post Internally round	Pack of 5	POM	TSPA 5	462030	В
For impressions on ground-down implant h that the impression can be poured with die s	-	·		ower limit of t	ne head, so



Double analogue	1 piece	Metal	IA4/IAU	462112	А
Castable abutment for large head Internally round	Pack of 5		РОВ	462086	В

CEMENTABLE ANGULATION ADAPTER (TI6AL4V)

These adapters are mounted on KOS **ROOT**[®] implants to compensate for the insertion direction. Plastic cements are preferably used. The implant head must be roughened beforehand. The protruding head parts are then removed. The impression is taken directly on the adapter.



TW2

Description	Code	REF	Price cat.
Adapter 15°	AA5 15°	462052	C
Adapter 25°	AA5 25°	462053	С

Description



INSERTION TOOLS





Length 23 mm Use with contra-angle	IT2W	900039	E
Use with ratchet	IT2 BCS	900030	E
Torque wrench 10-70 Ncm	TW2	425402	S

Code

REF

Price cat.

HEATLESS® DRILLS

Surgical steel, colour-coded, depth-coded and autoclaveable. The drill is marked with laser depth markings. Use between 3,000 and 5,000 rpm with good cooling and intermittent drill technique. Due to the extremely high cutting performance, you can work without pressure.



DRILLS MUST BE CHOOSEN DEPENDING ON BONE QUALITY AND IMPLANT DIAMETER.

Description	Colour	Max. working length	REF	Price cat.
DOS 1	yellow	17 mm	455311	D
DOS 2	black	17 mm	455312	D
DOS 3	red	17 mm	455313	D
DOS 4	blue	21 mm	455314	D
DOS 5	green	17 mm	455315	D
DOS 6	transparent	15 mm	455316	D

DOS 6 This drill is 2 mm shorter at the tip. It can therefore drill up to 2 mm deeper into hard bone than nominally indicated on the drill. Therefore, the conical bone cavity is only circularly extended in the crestal area without increasing the drilling depth.

PATHFINDER DRILLS

Conical 3-edge drill as initial drill, ideally suited for all crestal implant systems. The drill also passes between narrow cortical areas without pressure. For implants up to 20 mm length.

SUDX 1		
H		
-		
L	47 mm	I

Description	Colour	Max. working length	REF	Price cat.
BCDX 1	yellow	15 mm	900243	С
BCDX 2	black	15 mm	900244	С
BCDX 3	red	15 mm	900245	С

STARTER TRAY

Autoclaveable up to 134° C, not suitable for dry heat sterilizers. This surgical kit contains all drills and tools for first works with the system KOS **ROOT**[®]. Material: autoclaveable plastic.



on request

SURGERY

KOS ROOT[®] implants are used as compression screws. If possible, the hole should be created substantially thinner than the core diameter of the implant, since only in this way can good bone condensation be achieved. The minimum hole diameter depends on the bone density. For this reason, it is not possible to specify drill sequences that can be used favorably for all bone qualities. As a rule, it is necessary to drill much less into the soft maxilla (e.g. the DOS1 drill only can be used for KOS ROOT® implants with diameter 3.0-5.0) than into the well-mineralized mandible, which requires the use of a drilling sequence adjusted to the bone density. In very hard bone the implants should be inserted slighty deeper and then turned back 1/2 round.

1. Handling


Hold the implant by the holder and place the insertion tool on the implant head. The endosseous implant surface must not be touched. Pull out the implant with the plug and then twist off the plug with the needle holder at the predetermined breaking point.

KOS **ROOT**[®] implant with insertion tool **IT2W** (for angle piece) and **IT2 BCS** (manual)




Twisting off the bracket with the needle holder

3. Impression



- Attachment of the impression post
- Pressureless impression taking



Removal of the scoop from the implant. The impression post is located in the impression material. The impression can be sent to the laboratory.

2. Insertion using manual tools

Insert the implant by hand until it is firmly seated in the jaw.



Using the ratchet, torque ratchet or contraangle, screw the implant clockwise into the cavity. The endosseous part of the implant must be completely covered by bone. The polished implant neck is located in the mucosa. We recommend screwing the implant into the bone up to 1 mm deeper into the implant neck.

4. Laboratory procedures

Attachment of the impression post onto lab analogue.



IA4/IAU

The modeling is performed on the castable part POB (internally round; for bridges and bars).

Basal implants may only be used and operated by qualified persons with valid authorisation (para. 2 MedProdAnw Verordnung). We are certified according to DIN EN ISO 13485 and Annex II of Directive 93/42 EEC.

The product dimensions shown in this brochure may differ from reality for technical reasons.

KOS® is a registered trademark. Pat. Pend.

If implants are reprocessed, there is a risk of the development of infections, because no validated method for processing exists. Implants therefore may not be reprocessed.

COMPILATION AND EXPLANATION OF SYMBOLS ON THE PACKAGING:



STERILE R

Batch No.



Store tightly

keep closed

Non-sterile



Store in a dry

place





Do not use if packing

is damaged



Rx ONLY

Intended for use

by dentists or

surgeons only

Do not

resterilize



Manufacturer

Single use

product

1

Instruction

for use

Production

date



Catalogue number

Expiry date

(The products of this catalogue are CE marked (class I) and CE 1936 marked (class IIa and IIb) according to 93/42/EC Directive).

Commercial products that are not monitored by our notified body are declared as third-party products.

IHDEDENTAL*

Dr. Ihde Dental AG

Dorfplatz 11 CH - 8737 Gommiswald / SG Tel +41 (0)55 293 23 23 contact@implant.com www.implant.com

Dr. Ihde **Dental GmbH**

Erfurter Str. 19 D - 85386 Eching / Munich Tel +49 (0)89 319 761 0 Fax +49 (0)89 319 761 33 info@ihde-dental.de www.ihde-dental.de

