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### Hinweis

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## Two-part implant systems

### Note

Knowledge of the extensive specialized literature is prerequisite for the performance of implantological treatments. Surgical experience in the evaluation of the normal and pathological treatment is mandatory. The prosthetic planning must be performed prior to implantation - also if the surgical situation may cause a later alteration in the planning. These instructions alone are not sufficient for the immediate use of the Implant System.

### Description

The Implant System is a system for enossal dental implantations. It consists of implants, surgical and prosthetic instruments as well as abutments. The Implant System is suitable for single and two step implantation methods.

### Constituents

The individual active substances can be found on the respective package label or in the brochure concerning the use of the system.

### Indications

Upon adequate present bone material in respect of the bone quality, the bone width and the bone height:

- Single tooth replacement for spaces up to a width of incisors and premolars.
- Bridge replacements in dentate, partially dentate and non-dentate jaws
- Anchoring of prostheses via bridge and ball head systems
- Anchoring for orthodontic regulating systems

### Side effects

#### Temporary discomforts

Pain, feeling of tension, swelling, difficulties speaking and inflammation of the gingiva.

#### Long term disorders

- Oroantral and oronasal fistula
- Localized or systemic infections
- Crater shaped bone degeneration over the duration of use
- A horizontal and vertical bone degeneration is to be expected, also in proper surgical and prosthetic proceedings - as in all other dental implants. Type and extent of the bone loss is not foreseeable. Functionally caused and inflammatory bone degeneration can accumulate. If implants are placed in the vicinity of certain structures (nerves, maxillary sinus, neighboring teeth etc), a reversible or irreversible impairment (damage) of these structures can occur. The necessity of an antibiotic prophylaxis should be taken into consideration to prevent infections. Chronic pain has been reported in connection with dental implants. If an osseointegration does not occur, exfoliation of the implant must be expected.

### Contraindications

#### Absolute

- Patients whose systemic condition does not allow an implantation (reference is made to the studies of relevant literature for detailed understanding of this extensive contraindication)
- Prior to, during or after radiation therapy and / or malignant disorders
- Psychologic disorders, pain syndrome
- Uncompensated diabetes or other uncompensated systemic disorders
- Current or past therapy with bisphosphonates
- Disorders of the bone metabolism
- Bleeding tendency, immune deficiencies,
- Immediately loaded single tooth implantations
- Allergies to components of the implant system
- Titanium allergy
- Insufficient oral hygiene, lack of compliance
- Single tooth restorations on implants with a diameter 3.4 mm or less

#### Relative

- Diabetes, bruxism, allergies, pregnancy
- Alcohol and drug abuse, smoking
- Previously radiated bones, therapy with anticoagulants
- Disorders of the temporomandibular joint (on the other hand, numerous disorders of the temporomandibular joint can be treated with an implant therapy)
- Root remnants, periodontitis
- Limited available bone: In particular in the vertical dimension and / or in threatening vicinity of endangered structures (N. Metallis, maxillary sinus etc).
- If must always be considered if a basal implant (i.e. BOI) has less risks.

### Interactions

The interactions between titanium surfaces and oral hygiene products containing fluoride has been described in the literature. The patient is to be informed thereof. Due to the different rigidity of crestal implants and natural teeth, a greatly differing influence of the static is to be expected in combined works. Extensive implant prosthetic systems on crestal implants can impair the flexion and mobility of the cranial bones against one another. Pain conditions have been described in this connection.

### Use

#### Diagnosis / Explanation

The diagnosis consists of the anamnesis, the clinical examination, the radiological examination with small format x-ray and orthopantomogram as well as CT examinations. It is recommended to have a medical check-up performed by an experienced medical doctor.

Implantations are optional interventions, which require extensive explanation to the patient. Treatment alternatives and results of the implantations as well as results of a possible failure must be explained. In particular, the importance of oral hygiene and the costs of a following treatment must be explained. For type and extent of the explanation, reference is made for the respective case law.

### Approach

Under local anesthesia, preferably under flap formation, an implant bed is established with rotating or displacing instruments. For preparation of an adequate bone bed, the use of the drill and the compliance with a drilling sequence is recommended:

A corresponding additional distance to endangered structures must be complied with. For a detailed scheme concerning this aspect, please refer to our elaborate product brochures.

During the preparation of the bone cavity, sufficient cooling of the implant bed is to be ensured (i.e. with cooled (sterile) physiological saline solution).

Never drill above 500-800 rpm. Always use a sharp drill (max. 10 applications). Use intermittent reaming and drilling technique. The implant is removed from the sterile packaging immediately prior to introduction and stably inserted in the bone bed. It is mandatory that a good primary stability is achieved. Screw implants are screwed in with the help of a ratchet (RAT2), the torque wrench (TW2) or a screwdriver.

To minimize frictional warmth, it must be screwed in slowly and without great pressure. Depending on the available bone and the intended planning, the implants are inserted with the neck portion (polished shaft) extending from the bone or completely covered by the bone. If a removable dental prosthesis will be worn over the implants covered with mucous membrane during the healing phase, the implant must be inserted in the bone in order to avoid exposure during the healing phase as well to prevent chewing pressure on the implant. In the selection of the implant, the height and width as well as the quality of the present bone must be taken into consideration. It is recommended not to remove too much soft bone during the drilling, but rather to compress the bone to the side with displacing osteotomes.

Afterwards, the implant is either closed with the surgical screw on the mouth side (two step method), or a healing screw in adequate length (mucosa height) is screwed in (single step method). If indicated, an adequate wound closure and radiological control follow. During the healing phase (3-6 months are recommended depending on bone quality), the surgical situs is to be controlled in regular periods. In incomplete healing of the mucous membrane and / or premature mucous membrane penetration, it is recommended to replace the surgical screw with a healing cap. As long as the static of the total construction as well as the situation during the insertion allows, implants can also be loaded directly after the insertion. This is decided by the treating dentist. In case of doubt, a healing period of 3-5 months without loading should be complied with - depending on the quality of the bone. After completion of the healing phase, a radiological and clinical control is performed.

Abutments from the Implant System are available for the prosthetic treatment. The table below provides information concerning the fitting parts. For further information please refer to the respective brochure concerning the use of the system. Basal implants may only be handled and used by dentists / surgeons with valid authorisation issued by the manufacturer. This restriction includes all advices for patients, creating treatment plans, placement of implants, prosthetic treatments on implants, postoperative care (throughout the following years).

### Note

- After the implantation, the type of the used implant and the charge number must be noted in the patient file in writing.
- Implants may only be used in the period prior to expiration.
- Implants must be stored dry and sealed. The seal is to be opened immediately prior to introduction. All contact with foreign substances prior to insertion of the implant is to be avoided. The enossal part of the implant may not be touched.
- Upon accidental swallowing of implants, abutments or accessories, the remainder of the object is to be secured (x-ray or similar) and the necessary medical measures taken.
- The healing of the implant must be controlled regularly (i.e. radiological).
- To prevent bone loss, implants which are not osseointegrated or which are infected must be removed. The time of removal is determined by the dentist.

### Shelf life / Storage of the implants

The shelf life is 5 or 6 years after sterilization; see expiration date on the package. Store tightly closed and dry.

### Comment

Only originally sealed implants are sterile during the shelf life period. If implants are resterilized by the final consumer, no responsibility is assumed - independent of the method of sterilization. Sterile products are marked with STERILE. The expiration date is noted with the hour glass symbol. The LOT states the charge number. Implants may only be used once.

Dr. Ihde Dental AG reserves the right to change the design of the product and components or the packaging thereof, to adapt instructions of use as well as to newly agree on prices or delivery conditions. The liability is limited to the replacement of the defective product. Further claims of all types are excluded.

### Legend

	Read instructions
	Expiration date
	(Gamma) -sterilized
	Only use once
	Do not resterilize
	This product is only for sale to dentists, oral surgeons and maxillo-facial surgeons or on their behalf.
	Charge number
	Store in dry place
	Store tightly keep closed
	Temperature range from -5° C to 25° C
	Do not use if packing is damaged
	Manufacturer

CE 0029

996-2TEILIG-03

Revision issued 2016-01

**In case that implants would be reprocessed (cleaned, resterilized) infections could occur, because no validated procedures for reprocessing are available.**

**When processing resterilizable medical devices, the rules of the RKI-Guideline and the European standard EN ISO 17664 are relevant and must be obeyed. Guidelines regarding the processing of our instruments are available in the brochure accompanying our instruments. At the same time these guidelines are displayed on the website [www.implant.com](http://www.implant.com).**

**Keep this and all medicines out of the reach of children. To be used by dentists or surgeons only.**

**It is necessary that you carefully read the system application prospectus (SAP) for the implant system which you are planning to use. This prospectus contains important information about the application of the system's components.**

 **Dr. Ihde Dental AG**  
Dorfplatz 11  
CH-8737 Gommiswald/SG  
Tel. +41 (0)55 293 23 23  
Fax +41 (0)55 293 23 00  
contact@implant.com  
www.implant.com

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**Dr. InDe Dental AG**