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English: Instructions for use – Please read carefully

IHDENT® Chromat® Safeprint®

Indications for use: alginate for dental impressions.

CHROMAT: fast setting dust free chromatic alginate with color variations clearly indicating the critical phases of the impression taking: red during mixing, orange during working phase and positioning on the tray, yellow during the time in mouth while setting. Suitable in general clinical practice for the creation of study models. Aroma: Orange/Strawberry

Manufacturer Identification: Major Prodotti Dentari S.p.A. – Via Einaudi 23, 10024 Moncalieri (TO) - Italy

Dosing accessories: water beaker and powder spoon.

Procedure: Avoid moisture contamination: open the refill package and transfer its contents into an air tight container like MAJOR ALGINBOX immediately after opening. Before any use, shake appropriately the container to loosen the alginate powder.

Mixing ratio: 9,5 g powder : 20 ml water

Measuring: 1 (2) (3) (4) Spoon 9,5 g (19 g) (28,5 g) (38,0 g) powder: 20,0 ml (40,0 ml) (60,0 ml) (80,0 ml) water.

Directions: Dispense CHROMAT alginate powder from the container using the spoon provided and level the measure with a spatula. Place the powder into a mixing bowl and add the corresponding water. Two spoons of CHROMAT alginate and the corresponding amount of water are generally required for a complete impression. Mix the powder and water to a homogeneous consistency using a mixing spatula.

Refer to the manufacturer's instructions when using automatic mixers.

	CHROMAT
Mixing time:	35 sec.
Working time: *)	80 sec.
Minimum setting time in the mouth:	40 sec.
Setting time: *)	120 sec.

Times with *) apply from the start of mixing at a temperature of 23 °C/73 °F for both powder and water with deionized water. These times are reduced by higher temperatures and lengthened by lower temperatures. Variations in water hardness may also result in deviations from these times (if using tap water working and setting time are faster).

Model fabrication: Best results are achieved by rinsing the impression thoroughly after removal from the mouth, pouring immediately after cleaning for maximum accuracy. If longer storage is required, the impression should be kept in a firmly sealed plastic bag (relative humidity of 100 %).

This is essential for optimum results at a later pouring. The impression has a dimensional stability of 5 days. This means that the impression can be stored for 100 hours (in a firmly sealed plastic bag) without significant shrinking.

Disinfection: Impressions may be treated with standard disinfectant solutions. Please follow the respective manufacturer's instructions for disinfection and processing.

Product storage: CHROMAT alginates should be stored in a closed air tight container (like MAJOR ALGINBOX) in a dry area at maximum 25 °C/77 °F.

Shelf-life: As printed on packaging, our CHROMAT alginates have a shelf life of 6 years from manufacturing date.

Safety information for operators

The product is not subject to the classification and labeling requirements for hazardous mixtures as a medical device in its final state, destined for professional use (EU Regulation No. 1272/2008, Art. 1.5.d). The product is, therefore, also excluded from obligations of information (Material Safety Data Sheet, downstream safety information), as it is excluded from the application of Title IV EU Regulation No. 1907/2006 Art. 2.6.c.

This preparation contains dust that up to 5% may be composed of respirable crystalline silica (Cristobalite, CAS 144464-64-1); prolonged and massive inhalation of crystalline silica respirable dust may cause pulmonary fibrosis. In applying the following instructions, the quantity of dust released during daily dental practice should be such as to ensure that exposure levels below the risk thresholds are maintained. The material is, nevertheless, pretreated with binding agents in order to avoid the release of dust during measuring out and mixing processes.

Precautions

When measuring out the material, avoid dust formation and dispersion in the air. Do not inhale any dust that may be released. Close the containers immediately after use.

In the event of exposure.

When faced with a massive inhalation, blow your nose to facilitate respiration. If the dust gets in your eyes, do not rub or scratch them; rinse thoroughly with water or eye irrigating solution. The dust may cause temporary irritation of the respiratory system and skin. There are no other known significant dangers.

Occupational hygiene advice

Ensure that the measuring and mixing area is well-ventilated. Collect any material that may have been spilled with mechanical or suction devices. If mixed with water, the material may form slippery films.

Occupational exposure limits

The preparation does not contain powdered ingredients for which European Union exposure limits in the workplace exist. However, such limits do exist in some countries in Europe and around the world. For information, the following are some of the occupational exposure limits highlighted:

Exposure limits/ ingredients	Quantity	In some EU countries	In the USA (AC-GIH 2013)
General dust		TLV 10 mg/m ³ (inhalatable) TLV 3 mg/m ³ (respirable)	
Crystalline silica — cristobalite (in diatomaceous earth, soda ash flux calcined)	< 5.0% (respirable fraction)	0.05 mg/m ³ or above	TWA 0.025 mg/m ³
Zinc oxide	< 2.0%	2 mg/m ³ (Zn, inhaled) 0.1 mg/m ³ (Zn, respirable)	TWA 2 mg/m ³ (respirable) STEL 10 mg/m ³ (respirable)
Calcium sulfate dihydrate	< 20.0%	4 mg/m ³ (inhalatable) 1.5 mg/m ³ (respirable)	TWA 10.0 mg/m ³ (inhalatable)
Talc	< 0.5%	0.5 mg/m ³ or above	2.0 mg/m ³ (inhalatable)

Note

The product is not subject to the classification and labeling requirements for hazardous mixtures as a medical device in its final state, destined for professional use (EU Regulation No. 1272/2008, Art. 1.5.d). The product is, therefore, also excluded from obligations of information (Material Safety Data Sheet, downstream safety information), as it is excluded from the application of Title IV EU Regulation No. 1907/2006 Art. 2.6.c. Therefore the Material Safety Data Sheet is not available in accordance with Art. 31 or 32 of the stipulated EU Regulation No. 1907/2006.

9999-7600

Revision issued 2016-06

Legends

- Read instructions
- Caution, consult accompanying documents
- LOT Charge number
- Keep in a dry place
- Der grüne Punkt
- Manufacturer

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IHDENT® is a registered trademark.
Dr. Ihde Dental AG
Dr. Ihde Dental GmbH