Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Calciumhydroxid liner

Calcium hydroxide is a radio-opaque, low viscosity liner with an aqueous base. Calcium hydroxide is administered in application syringes or tabs.

Significant components

Calcium hydroxide

Areas of use

- Indirect pulpa capping for the treatment of caries profunda.
- direct pulpa capping.
- Underfilling of cavities in order to prevent the exposition of acidic mediums in the use of cements (stone cement).
- Temporary fillings of the root canal if the growth of the root and the apex formation is not yet completed.

Use

Clean and dry the prepared cavity. If necessary, the application needle can be individually formed to a slight arc in order to ease the application.

Note

The first drop of the liner is to be disposed of. Air encapsulated during filling can otherwise cause an imposition of the material. This also applies after each exchange of the application needle.

Apply calcium hydroxide exactly in a thin layer. The applied layer is dried with a water and oil free air flow. The substance may not be applied to the edges of the cavity. Paste remains are therefore to be completely removed. After completion of the application and drying, continue with the processing.

Cover fillings can be applied as composites, compomeres, glass ionomer cements or amalgam.

For filling the root, the root canal is prepared and cleansed (e.g. with 3 % hydrogen peroxide or sodium hypochlorite solution). Afterwards the caries profund calcium hydroxide is filled into the canal with a suitable instrument. Avoid over filling (apical) and encapsulation of air. It is recommended to check the root canal filling with x-ray imaging.

Side effects

If the root canal is over filled, unpleasant side effects can occur such as sensitivity to pressure and/or swelling. These normally disappear within 24 - 48 hours. The necessity of an intervention should be carefully considered in such cases.

Contraindications or interactions

not known

Precautionary statements

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P280 Wear protective gloves / protective clothing / eye protection / face protection.

Store out of reach of children. Only for use by the dentist. Also follow the warning in the safety data sheets. These can be sent upon request.

Storage

If the needle should remain on the syringe after use, the needle must be inserted through the membrane into the glass container in order to avoid drying of the calcium hydroxide. The needle can alternatively also be removed and the provided seal screwed on the tip. A new needle is to be attached to the syringe for the next use. **Do not use hardened paste, therefor do not use tabs twice.**

Store dry, not below 8 °C and max. 25 °C. See barcode on the package for charge number and expiration date.

Package size

Application syringe with 2 ml content (approx. 3,5 g)

Accessories: 3 application needles,

storage container with septum against drying

REF 8768

Replacement tips, bag of 3 pieces

REF 4058

Legends

Read instructions

Expiration date

LOT Charge number

For single use only

Keep in a dry place

Store tightly closed keep

Temperature range from 8 °C to 25 °C

Manufacturer

CE 0029

Danger

H315 Causes skin irritation.

H318 Causes serious eye damage.

H335 May cause respiratory irritation.

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IHDENT® is a registered trademark.
**Dr. Ihde Dental AG**
Dr. Ihde Dental GmbH