Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® P&B®

IHDENT® P&B® is a light curing priming and bonding system that is easy to use and which does not require mixing prior to application. It was developed for a strong binding of composites and amalgam to tooth enamel and dentin as well as composites of non-precious and precious metals.

It's strong adhesion to tooth enamel or dentin occurs in a similar manner as with glass ionomer cements. A good, durable adhesive strength and good biocompatibility are achieved with the methacrylic carbon acid ester. IHDENT® P&B® therefore differs significantly from other adhesive systems, with which the adhesion is achieved with phosphate groups, which are however sensitive to hydrolysis and the adhesive strength can therefore become weaker over time.

IHDENT® P&B® is compatible with all light curing composite materials which are currently available on the market.

Constituents

Mixture of mono and difunctional methacrylates, polyacrylic acids, acetone, photo initiators, accelerators and stabilizers.

Side effects

Avoid contact with the skin, the mucous membranes and the eyes since allergies cannot be ruled out.

Contraindications

Avoid direct contact with the pulpa. IHDENT® P&B® is not to be used in persons who are known to react allergic to methacrylate.

Interactions

Materials containing eugenol are not to be used in connection with IHDENT® P&B® since these could influence the curing of the material.

Therapeutic indications

Adhesive of light curing composites and amalgam on tooth enamel and dentin as well as compositions of non-precious and precious metals.

Use

To assure isolation against moisture, the use of a cofferdam is recommended.

1. Preparation of the cavity

Prior to the preparation, the tooth is to be cleansed of residues, concrements etc. The cavity is to be prepared so that as little possible of the healthy tooth material is lost. The cavity margin is to be a little chamfered (0.5 - 1.0 mm) on the tooth enamel so that the adhesive surface of the tooth enamel and thereby the adhesive strength can be increased.

2. Protection of the pulpa

In deep drillings, the cavity floor close to the pulpa must be covered with a thin layer of curing calcium-hydroxide underfilling material.

3. Enamel and dentin conditioning

The total etching technique is recommended: ETCHING GEL (REF 5138) is applied to the enamel and dentin surfaces, starting with the enamel margins and allowed to work for 20 seconds. In deciduous teeth or highly fluoridated teeth, an etching time of 60 seconds is recommended. Afterwards rinse with water for 20 seconds and dry with water and oil free air. The dentin should not be overdried, since a slightly damp dentin surface is important for the function of IHDENT® P&B®. The etched enamel margin must have a chalk white appearance.

Precautionary measures during etching: It is important that tooth enamel already etched is not contaminated. Should a contamination occur, the etching must be repeated as described above. Avoid contact of the etching gel with the skin, mucous membranes and eyes. Should accidental contact occur, rinse thoroughly with water.

4. Application of IHDENT® P&B®

A sufficient amount of IHDENT® P&B® is applied to the dentin and enamel surfaces with a brush and the primer is intensively brushed in for 30 seconds. Note is to be taken that the dentin surfaces to be impregnated are kept damp with the primer over the entire prescribed time.

Attention is to be taken that the material is distributed homogenous and over the entire surface. Remove possible excess material. All volatile components are removed afterwards through soft blowing with pressurized air (approx. 15 seconds) and the adhesive thereby distributed to an even layer. Do not overdry the dentin.

Note: Do not spill IHDENT® P&B®! Close the bottle immediately after the removal of IHDENT® P&B®. If it is not used immediately, IHDENT® P&B® which has already been removed from the bottle is to be stored under low lighting in order to avoid premature polymerization from light.

5. Curing of IHDENT® P&B®

The IHDENT® P&B® coating is cured with 20 seconds of constant radiation with a suitable dental lamp.

Note: IHDENT® P&B® does not cure by itself. It is important that the primed dentin surface and the etched enamel are kept dry and clean before the next application is performed.

6. Application of a second layer of IHDENT® P&B®

The application of a second layer of IHDENT® P&B® is recommended. This is done in the same manner as the first application as described above. After curing of this second layer with a dental lamp, a light curing composite can be applied.

7. Laying the filling, curing and finishing

Please refer to the instructions for the light curing filling material used.

Warning: IHDENT® P&B® contains acetone. Acetone is highly flammable. Keep away from sources of ignition. Do not smoke. Do not inhale fumes. Take measures against electrostatic charging.

Precautionary statements

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Keep out of reach of children. Only to be used by the dentist.

Storage instructions

Do not store over 23 °C (74° F)! Store material in a well ventilated location. Do not use after the expiration date! See barcode on the package or bottle label for charge no. and expiration date.

Package sizes

10 ml IHDENT® P&B®	REF 5120
3 ml IHDENT® P&B®	REF 5121

Legends

	Read instructions
	Expiration date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperatures up to 23 °C
	Manufacturer

CE 0029

Danger

H225 Highly flammable liquid and vapour.

H319 Causes serious eye irritation.

H336 May cause drowsiness or dizziness.

9999-5120-03

Revision issued

2015-07

Dr. Ihde Dental AG
Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH
Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT® and P&B® are registered trademarks.

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