Deutsch:

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English:

Instructions for use – Please read carefully

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IHDENT® Safeprint® IM is a syringeable, monophase precision impression material for partial or full impressions for crowns, bridges, inlays, partials or full dentures as well as dental repairs. IHDENT® Safeprint® IM is also suitable as an impression material for reline impressions, functional impressions, denture repair impressions and for the copper-ring technique. IHDENT® Safeprint® IM can also be used in conjunction with IHDENT® Safeprint® P in the two-step and one-step technique.

IHDENT® Safeprint® IM has a medium to high viscosity and is based on addition-curing silicones. It is supplied in automatic-mixing cartridges and can be applied directly without bubbles to the tooth or impression tray. IHDENT® Safeprint® IM has hydrophilic properties, which enhance the wettability and reproduction of detail, particularly in the subgingival regions. IHDENT® Safeprint® IM is thixotropic and is therefore positionally stable in the impression tray, but is flowable when pressure is applied. IHDENT® Safeprint® IM meets the requirements of DIN 13913 1-b and ISO 4823 Cat. A, Type 2 (medium viscosity).

## Preparing the dispensing syringe

Remove the seal of the AUTOMIX double cartridge and replace it with an AUTOMIX mixing tip. Then fit the AUTOMIX cartridge to the dispensing gun and the dispensing gun is now ready for application. Additional intraoral tips are intended for direct, fine application to the tooth.

## Taking the impression

Load the automatically mixed material onto the tray or apply it directly to the tooth. The impression tray should be placed in the patient's mouth after no longer than 2 minutes and 15 seconds and should remain in situ for at least 3.5 minutes.

## Pouring the impression

The impression should be stored at normal room temperature (max. 25 °C). The impression can be poured 30 minutes after removal from the mouth. There are no further time restrictions.

## Electroforming

All impressions can be copper or silver electroformed.

## Disinfection

Impressions can be disinfected using commercially available disinfectant solutions. (We recommend IHDENT® Dentatrend® Spray disinfection REF D0102C).

## Additional information

Do not use with standard silicone impression materials (condensation reaction type).

## Processing

IHDENT® Safeprint® IM in **cartridges**:

|  |                       |
|--|-----------------------|
| Mixing time                            | n/a                   |
| Total working time (incl. mixing time) | <b>2 min. 15 sec.</b> |
| Setting time (incl. mixing time)       | <b>4 min.</b>         |
| Recommended time in situ               | <b>3 min. 30 sec.</b> |
| Dimensional change after 24 hrs        | <b>0.1 %</b>          |
| Permanent set                          | <b>0.4 %</b>          |
| Compression set                        | <b>3.0 %</b>          |
| Shore A hardness                       | <b>approx. 62</b>     |

## Side effects

None known.

## Interactions

Do not use latex gloves, as they may contaminate the material and impair the setting reaction of the polyvinyls.

## Contraindications

Do not use with condensation-curing silicone impression materials!

## Storage

IHDENT® Safeprint® IM has a minimum shelf life of 3 years, provided it is stored at 23 °C. Do not store over 25 °C. Long-term storage of the sealed pack shortens the working and setting times. Refer to the barcode on the packaging and cartridge label for the Batch No. and Expiry date.

**Keep out of reach of children. Intended for use by dentists only.**

## Package size

2 cartridges à 50 ml, Base- and catalypaster,  
12 mixing cannulas, pink  
for safety mixing gun (1:1) **REF 7755**

## Legends

|  |                           |
|--|---------------------------|
|  | Read instructions         |
|  | Expiration date           |
|  | LOT Charge number         |
|  | Keep in a dry place       |
|  | Store tightly closed keep |
|  | Temperatures up to 25 °C  |
|  | Manufacturer              |

**DIN EN ISO 4823:2007**

9999-7755-03

Revision issued 2015-08

**Dr. Ihde Dental AG**  
Dorfplatz 11  
CH-8737 Gommiswald/SG  
Tel. +41 (0)55 293 23 23  
Fax +41 (0)55 293 23 00  
contact@implant.com  
www.implant.com

## Distribution by

**Dr. Ihde Dental GmbH**  
Erfurter Str. 19  
D-85386 Eching/München  
Tel. +49 (0)89 319 761-0  
Fax +49 (0)89 319 761-33  
info@ihde-dental.de  
www.ihde-dental.de

IHDENT® and Safeprint® are registered trademarks.

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