Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Safeprint® P / LV / MV

Addition-cured silicone precision impression material (polyvinyl siloxane)

IHDENT® Safeprint® is available in 3 consistencies to facilitate use with any standard dental impression technique. The optimal mixing ratio of 1:1 (volume or weight) ensures reliable, accurate dispensing. The different pigmentation of the base and catalyst materials guarantees a reliable check of the mixture.

IHDENT® Safeprint® P is a putty impression material and is used in conjunction with Safeprint® LV for the one-step putty-wash impression technique; it complies with EN ISO 4823:2007 Type 0 (putty). Safeprint® P has low hydrophilic properties, is easily mixed and does not stick.

IHDENT® Safeprint® LV is a low viscosity wash material (syringe material) for the one-step putty-wash impression technique (e.g. with Safeprint® P) and complies with DIN EN ISO 4823:2007 Type 3 (low viscosity).

IHDENT® Safeprint® MV is a syringeable Monophase precision impression material for crowns, bridges, inlays, partial and full dentures as well as dental repairs. IHDENT® Safeprint® MV is also suitable as an impression material for reline impressions, functional impressions, denture repair impressions and for the copper ring technique as well as a syringe material for two-step direct impressions (e.g. with IHDENT® Safeprint® P). IHDENT® Safeprint® MV complies with DIN EN ISO 4823:2007 Type 3 (low viscosity).

Application

IHDENT® Safeprint® P: Perforated or non-perforated metal or plastic impression trays can be used. Trays should be coated with a thin film of silicone adhesive to improve the bond of IHDENT® Safeprint® P to the trays. Allow the adhesive to dry for 3 minutes. Remove equal amounts of impression material with the measuring scoops (white for the catalyst, black for the base) and knead for a minimum of 30 seconds to attain a uniform colour. Load the mixed material onto the prepared tray and place into the patient's mouth no longer than 1:30 minutes (working time at 23 °C) after the beginning of mixing. Allow to cure intraorally for a minimum of 3 minutes.

IHDENT® Safeprint® LV and MV

Preparing the mixing gun

Remove the seal of the safety cartridge and replace it with a mixing tip. Place the cartridge into the mixing gun. Additional intraoral tips are intended for direct, fine application to the tooth.

Taking the impression

Load the automatically mixed material onto the impression tray or apply it directly to the tooth. The impression tray should be placed in the patient's mouth after no longer than 2 minutes (working time at 23 °C) and should remain in situ for at least 3 min.

Note

The used mixing tip should remain fitted as a seal until the next application. **Do not replace the original seal on the safety cartridge**, as this may cause cross contamination and unintentional partial curing cannot be excluded.

Application	Safeprint® P	Safeprint® LV	Safeprint® MV
Consistency	putty	low viscosity	medium viscosity
Mixing time*	at least 30 sec.	at least 30 sec.	at least 30 sec.
Total working time*	1.75 min.	2.25 min.	2.25 min.
Intraoral setting time*	3 min.	3 min.	3 min.
Setting time*	5 min.	5 min.	5 min.
Dispensing base/catalyst	1:1	1:1	1:1
Dispensing aids	measuring scoop white and black	none	none

* At a room temperature of 23 °C. Higher temperatures shorten the prescribed setting time and lower temperatures extend the prescribed setting time.

Pouring the impression

The impression should be stored at normal room temperature (max. 25 °C). Allow a minimum of 30 minutes after removal from the mouth before pouring the impression. There are no further time restrictions.

Electroforming

All impressions can be copper or silver electroformed.

Disinfection

Impressions can be disinfected using standard disinfectant solutions.

Storage / Shelf life

Do not use the material after the expiry date. Refer to the adhesive label on the folding carton for the Batch No. and Expiry date.

Note

Do not interchange measuring scoops (IHDENT® Safeprint® P)! Do not use with standard silicone impression materials (condensation-cured)! Do not use latex gloves, as they may contaminate the material and impair the setting reaction of the polyvinyls.

Keep out of reach of children. Intended for use by dentists only.

Package sizes

IHDENT® Safeprint® P

1 container base paste (green), 1 container catalyst paste (white),
ever 450 ml = ever ca. 710 g **REF 7710**

IHDENT® Safeprint® LV

2 cartridges à 50 ml, 12 mixing cannulas, yellow **REF 7725**

IHDENT® Safeprint® MV

2 cartridges à 50 ml, 12 mixing cannulas, yellow **REF 7735**

Accessories IHDENT® Safeprint® LV & MV

Safety mixing gun (1:1) **REF 7746**

24 Mixing cannulas, yellow **REF 7767**

50 Intraoral tips, yellow **REF 7770**

Legends

-  Read instructions
-  Expiration date
-  LOT Charge number
-  Keep in a dry place
-  Store tightly closed keep
-  Temperatures up to 25 °C
-  Manufacturer

DIN EN ISO 4823:2007

9999-7710-03

Revision issued 2015-08

Dr. Ihde Dental AG
Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH
Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT® and Safeprint® are registered trademarks.

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